# PMSX SELF EXPANDING

- FR** Stent vasculaire périphérique auto-expansible
- EN** Peripheral vascular self expanding stent
- DE** Peripherer vaskulärer selbstexpandierender stent
- ES** Stent vascular periférico autoexpandible
- IT** Stent autoespandibile vascolare periferico
- PT** Stent auto-expansível vascular periférico
- NL** Perifeer vasculair zelfexpanderend stent
- PL** Stent samorozprężalny do naczyń obwodowych
- DA** Perifer vaskulær selvudvidende stent
- SV** Perifer självutvidgande stent
- NO** Selvekspanderende stent til perifere kar
- FI** Itselaajeneva ääreisverenkierron stentti
- EL** Περιφερικό αγγειακό αυτοεκτεινόμενο στεντ
- TR** Kendiliğinden genişleyen periferel vasküler stent
- RO** Stent autoexpandabil vascular periferic
- HU** Perifériás vaszkuláris öntáguló sztent rendszer
- RU** Периферийный васкулярный саморасширяющийся стент
- CS** Periferní vaskulární samoroztažitelný stent
- BG** Периферен саморазширяващ се стент
- AR** دعامة إستنت ذاتي التمدد للأوعية الدموية المحيطة

**PEROUSE**  
MEDICAL

**STERILE**

FR : Stérilisé avec du gaz d'oxyde d'éthylène  
EN : Sterilized with ethylene oxide gas  
DE : Mit Ethylenoxidgas sterilisiert  
ES : Esterilizado con óxido de etileno  
IT : Sterilizzato ad ossido di etilene  
PT : Esterilizado com gás óxido de etileno  
NL : Gesteriliseerd met ethyleenoxide  
PL : Poddano sterylizacji gazowej tlenkiem etylenu  
DA : Ethylenoxid gas sterilisation  
SV : Steriliserat med etylenoxid  
NO : Sterilisert med etylenoksidgass  
FI : Steriloitu etyleeni-oksidaasulla  
EL : Αποστειρωμένο με αέριο αιθυλενοξείδιο  
TR : Etilen oksit gazı ile sterilize edilmiştir  
RO : Sterilizat cu gaz oxid de etilenă  
HU : Etilén-oxid gázzal sterilizálva  
RU : Стерилизовано при помощи оксида этилена.  
CS : Sterilizováno plynným etylenoxidem  
BG : Стерилизирано с газообразен етиленов оксид  
AR : يتم تعقيمه باستخدام غاز أكسيد الإيثيلين

FR: Attention  
EN: Caution  
DE: Achtung  
ES: Precaución  
IT: Attenzione  
PT: Atenção!  
NL: Let op  
PL: Przestroga  
DA: Advarsel  
SV: Fara  
NO: Vær forsiktig  
FI: Varoitus  
EL: Προσοχή  
TR: Dikkat  
RO: Atenție  
HU: Figyelmeztetés  
RU: Осторожно  
CS: Pozor  
BG: Внимание  
AR: تحذير

FR: Usage unique  
EN: Single use  
DE: Zum einmaligen Gebrauch  
ES: De un solo uso  
IT: Monouso  
PT: Utilização única  
NL: Eenmalig gebruik  
PL: Do jednorazowego użytku  
DA: Engangsbrug  
SV: Engangsbruk  
NO: Engangsbruk  
FI: Kertakäyttö  
EL: Μιας χρήσης  
TR: Tek kullanımlık  
RO: De unică folosință  
HU: Egyszeri használatra  
RU: Одноразового применения  
CS: Najedno použití  
BG: За еднократна употреба  
AR: يُستخدم مرة واحدة

FR: Ne pas restériliser  
EN: Do not reesterilize  
DE: Nicht erneut sterilisieren  
ES: No reesterilizar  
IT: Non risterilizzare  
PT: Não voltar a esterilizar  
NL: Niet opnieuw steriliseren  
PL: Nie sterylizować ponownie  
DA: Må ikke gensteriliseres  
SV: Får inte omsteriliseras  
NO: Må ikke reesteriliseres  
FI: Ei saa steriloida uudestaan  
EL: Μην επαναποστειρώνετε  
TR: Tekrar sterilize etmeyin  
RO: Nu sterilizați a doua oară  
HU: Ne sterilizálja újra!  
RU: Не стерилизовать повторно  
CS: Nesterilizujte opakovaně  
BG: Да не се стерилизира повторно  
AR: عدم إعادة تعقيمه

FR: Conserver à l'abri de la lumière du soleil  
EN: Keep away from sunlight  
DE: Unter Lichtabschluss aufbewahren  
ES: Conservar en un lugar protegido de la luz solar  
IT: Conservare al riparo dalla luce del sole  
PT: Manter abrigado da luz solar  
NL: Beschermen tegen het zonlicht  
PL: Przechowywać z daleka od słońca  
DA: Beskyttes mod sollys  
SV: Skyddas från solljus  
NO: Lagres beskyttet fra sollys  
FI: Säilytettävä suojassa auringonvalolta  
EL: Να φυλλάσσεται μακριά από το ηλιακό φως  
TR: Güneş ışığından uzak tutunuz  
RO: A se păstra la adăpost de lumina soarelui  
HU: Napfénytől távol tartandó  
RU: Хранить в месте, защищенном от солнечных лучей  
CS: Zabraňte přístupu slunečního světla  
BG: Да се пази от пряка слънчева светлина  
AR: الإطلاع على نشرة التعليمات

FR: Craint l'humidité  
EN: Protect from humidity  
DE: Vor Feuchtigkeit schützen  
ES: Proteger de la humedad  
IT: Teme l'umidità  
PT: Manter a humidade  
NL: Kan niet tegen vocht  
PL: Trzymać z daleka od wilgoci  
DA: Beskyttes mod fugt  
SV: Skyddas från fukt  
NO: Beskyttes mot fuktighet  
FI: Suojattava kosteudelta  
EL: Να προστατεύεται από την υγρασία  
TR: Nemi koruyunuz  
RO: A se proteja de umezeală  
HU: Száraz helyen tartandó  
RU: Поддерживайте влажность  
CS: Udržujte vlhkost  
BG: Да се пази от влага  
AR: يُحفظ بعيداً عن الرطوبة

FR: Consulter les précautions d'emploi  
EN: See Instructions For Use  
DE: Anleitung lesen  
ES: Consulte las instrucciones de uso  
IT: Leggere attentamente le istruzioni per l'uso  
PT: Consulte as instruções de utilização  
NL: Lees de handleiding  
PL: Zobacz instrukcję obsługi  
DA: Se brugsanvisning  
SV: Se bruksanvisningen  
NO: Les bruksanvisningen  
FI: Lue käyttöohje  
EL: Βλ. οδηγίες χρήσης  
TR: Kullanma kılavuzuna bakınız  
RO: Consultați instrucțiunile de utilizare  
HU: Lásd a használati útmutatót  
RU: См. инструкцию по использованию  
CS: Viz návod k použití  
BG: Вижте указанията  
AR: الإطلاع على نشرة التعليمات

FR: À utiliser avant la date  
EN: Use by  
DE: Verfallsdatum  
ES: Utilizar antes de  
IT: Utilizzare entro  
PT: Utilizar antes de  
NL: Te gebruiken vóór  
PL: Do użycia przed  
DA: Anvendes inden  
SV: Används före  
NO: Brukes før  
FI: Käytettävä ennen  
EL: Χρήση έως  
TR: Son kullanma tarihi  
RO: A se utiliza înainte de  
HU: Minőségét megőrzi:  
RU: Использовать до  
CS: Použitelné do  
BG: Да се използва преди  
AR: تاريخ إنتهاء الصلاحية

FR: Ne pas utiliser si l'emballage est endommagé  
EN: Do not use if the packaging is damaged  
DE: Nicht verwenden, wenn die Verpackung beschädigt ist  
ES: No utilizar si el embalaje está deteriorado  
IT: Non utilizzare se la confezione è danneggiata  
PT: Não utilizar se a embalagem se encontrar danificada  
NL: Niet gebruiken als de verpakking beschadigd is  
PL: Nie używać jeśli opakowanie jest uszkodzone  
DA: Må ikke anvendes, hvis emballagen er beskadiget  
SV: Får inte användas om förpackningen är skadad  
NO: Skal ikke brukes hvis emballasjen er skadet  
FI: Ei saa käyttää, jos pakkaus on vioittunut  
EL: Να μην χρησιμοποιείται αν η συσκευασία έχει υποστεί ζημιά  
TR: Ambalaj hasar görmüşse kullanmayınız  
RO: A nu se utiliza dacă ambalajul este deteriorat  
HU: Ne használja fel, ha a csomagolás sérült  
RU: Не использовать в случае повреждения упаковки  
CS: Neroužívejte, pokud je porušen obal  
BG: Да не се употребява, ако е нарушена целостта на опаковката.  
AR: لا يُستخدم في حالة وجود أي تلف بالعبوة

LOT

FR: Numéro de lot  
EN: Lot number  
DE: Chargennummer  
ES: Número de lote  
IT: Numero di lotto  
PT: Número de lote  
NL: Lotnummer  
PL: Numer partii produktu  
DA: Batch-kode  
SV: Lotnummer  
NO: Lotnummer  
FI: Eränumero  
EL: Κωδικός παρτίδας  
TR: Parti kodu  
RO: Număr lot  
HU: Adagszám  
RU: Номер партии  
CS: Číslo šarže  
BG: Партида номер  
AR: رقم التتبع

FR: Fabricant  
EN: Manufacturer  
DE: Hersteller  
ES: Fabricante  
IT: Produttore  
PT: Fabricante  
NL: Fabrikant  
PL: Producent  
DA: Fabrikant  
SV: Tillverkare  
NO: Manufacturer/Produsent  
FI: Valmistaja  
EL: Κατασκευαστής  
TR: Üretici  
RO: Producător  
HU: Gyártó  
RU: Производитель  
CS: Výrobce  
BG: Производител  
AR: الصانع


Figure 1 – Système de stent

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
Figure 1 – Stent system

**STERILE:** Sterilized with ethylene oxide. Do not use if the sterile package is open or damaged.

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions; otherwise, complications may occur.

Use the product before expiration of the “Best-Before-Date” mentioned on the packing.

## I. DESCRIPTION

**Please use the product illustrations above to guide you through the device description.**

The Peripheral Vascular Self Expanding Stent System is designed to deliver a self-expanding stent to the peripheral vasculature via a sheathed delivery device. The Peripheral Vascular Self Expanding Stent System is comprised of the following:

- An implantable self-expanding nickel-titanium alloy (nitinol) stent, as shown in Figure 2. The stent is flexible, a fine tubular mesh prosthesis, with a helical design, which achieves its unconstrained diameter upon deployment into the target vessel. Upon deployment, the stent imparts an outward radial force on the luminal surface of the vessel to establish patency.


Figure 2 – Nitinol Stent

- A stent system, as shown in Figure 1, comprised of an inner tubing assembly that contains the guide wire lumen, a stent delivery sheath and a system stability sheath, which are linked together by means of a handle. The guidewire lumen terminates distally in an atraumatic delivery device catheter tip and originates proximally in a luer hub designed to accept a compatible guide wire.
- The self-expanding stent is constrained in the space between the guide wire lumen and stent delivery sheath. Unintended stent movement during sheath retraction is restricted by the delivery device. The stent delivery sheath has a radiopaque zone at its distal end. The tip is radiopaque, too. Proximal and distal to the loaded stents are radiopaque markers. Prior to deployment the shipping lock (Figure 3) must be removed and discarded.


Figure 3 – Handling of shipping lock

## CONTENT

- 1 Stent system
- 1 Instruction for use

## II. RECOMMENDED STORAGE

1. Keep in a cool, dark and dry place.
2. Use immediately after the sterile package is opened.
3. Refer to the symbol legend at the end of this document.

## III. DISPOSAL INSTRUCTIONS

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

## IV. INDICATIONS

The Peripheral Vascular Self Expanding Stent System is intended for primary stenting of de-novo or restenotic lesions of the peripheral arteries.

## V. CONTRAINDICATIONS

Generally, contraindications to PTA are also contraindications for stent placement. Contraindications include, but are not limited to:

- Patients with highly calcified lesions resistant to PTA.
- Patients with a target lesion with a large amount of adjacent acute or subacute thrombus.
- Patients with Hyper-coagulopathy.
- Patients with stenosis or occlusion where lesion crossing with guide wire is not possible.
- Patients with fresh, soft thrombotic or embolic material.
- Patients with asymptomatic stenosis (Fontaine-Stage I).

## VI. WARNINGS AND PRECAUTIONS

1. The Peripheral Vascular Self Expanding Stent System is supplied sterile and is intended for single use only. DO NOT resterilize and/or reuse the device.
2. DO NOT use if pouch is opened or damaged.
3. DO NOT use the stent system after the end of the month indicated by the “Expires” date specified on the package.
4. Persons with allergic reactions to nickel titanium (nitinol) alloy may suffer an allergic response to this implant.
5. DO NOT use with ETIODOL™ or Lipiodol contrast media.
6. DO NOT expose the stent system to organic solvents (e.g. alcohol).
7. The stent is not designed for repositioning or recapturing.
8. Stenting across a major branch could cause difficulties during future diagnostic or therapeutic procedures.
9. This medical device may be used only by physicians skilled in percutaneous transluminal angioplasty.

10. Prior to use, the packaging and product should be inspected for signs of damage. DO NOT use if packaging is damaged.
11. DO NOT use agents containing organic solvents or oleaginous contrast media. Contact with these agents may lead to damage of the product.
12. The Peripheral Vascular Self Expanding Stent System should be handled carefully. Prior to use inspect the Peripheral Vascular Self Expanding Stent System carefully for bends, kinks, or other damage. DO NOT use a damaged Peripheral Vascular Self Expanding Stent System.
13. Confirm the compatibility of the product diameter and length with the introducer sheath and guide wire prior to use.
14. Since serious complications might arise when using this Self Expanding Stent System, operation should be done in a medical institution where emergency procedures can be executed.
15. This Peripheral Vascular Self Expanding Stent System can only be inserted with the use of a guide wire (Insertion of this device alone may lead to damage to the vascular wall or perforation of vessels.)
16. The guide wire must be completely advanced to reach towards the end of the vessel containing the lesion to be treated. [If guide wire is not completely advanced to the end, it may loosen from the guide wire lumen and it may even lead to damage to the vascular wall or perforation of vessels.]
17. If the guide wire is displaced from the guide wire lumen during operation, remove the device and reinsert the guide wire.
18. If abnormal or strong resistance is experienced during the operation the cause for such abnormality or resistance should be verified and appropriate measures should be performed before proceeding [If such abnormality or resistance is ignored and excessive force is applied, it may lead to damage of the vessels or to the device shaft breaking and remaining inside the body.]
19. During usage the self-expanding stent system should be replaced with a new system if there are any visible kinks or breaks in the device shaft that are found during handling or implantation.
20. If a great resistance is encountered during insertion, movement, or pulling out of the Peripheral Vascular Self Expanding Stent System it should be verified with fluoroscopy that the guide wire is not tangled around the distal end of the device and that the distal end of the devices is not caught on the stent.
21. This device should not be unreasonably inserted into or swiftly pulled out from lesions with highly tortuous vessels, bifurcated lesions, or calcified lesions. [The shaft at the distal end may then kink or be damaged, leading to vascular damage.]
22. Precaution should be taken to prevent any damage to the device by other equipment (such as scalpels, blades or scissors). DO NOT use a damaged device.
23. During the usage of this device, the temperature, blood pressure, pulse, and respiration of patients should be monitored. In case of any abnormality, the procedure should be stopped or appropriate measures taken based on the physician’s judgment.
24. Refer to the instructions supplied with any interventional devices to be used in conjunction with the system for their intended uses, contraindications, warnings, precautions, and instructions for use.
25. Appropriate drug therapy (anticoagulant, vasodilator, etc.) should be administered to the patient according to standard protocols for percutaneous interventions before insertion of the Peripheral Vascular Self Expanding Stent System.
26. The delivery device is not designed for use with power injection systems.
27. Recrossing a partially or fully deployed stent with adjunct devices must be performed with caution.
28. Prior to stent deployment, remove slack (see marked area in figure 4) from the stent system outside the patient.
29. DO NOT attempt to break, damage, or disrupt the stent after placement.
30. As under PTA normally bigger and/or longer sheaths are used, which remain normally for a longer time than during a standard diagnostic catheter, vascular complications can be caused (examination hematoma, arterial fistula, aneurysms spuria and arterial thromboses).
31. When retraction of the delivery device is not possible image with fluoroscopy to see if the device is trapped in any way. If there is no appearance of entrapment in any way slowly rotate the device handle 10 degrees in alternating direction while applying slight traction. If the device cannot be retrieved and the outer sheath is retracted after deployment carefully recover the stent pod area as this may allow the device to be removed easily through the vessel or guide catheter using the same or a similar technique as described above.

## VII. RECOMMENDED DRUG REGIMEN

Patients should receive adequate anti-platelet and anticoagulation therapy as prescribed by their physicians.

## VIII. POSSIBLE COMPLICATIONS AND ADVERSE REACTIONS

Potential complications associated with the use of peripheral stents may include, but are not limited to:

- Allergic/anaphylactoid reaction
- Amputation
- Aneurysm
- Angina/coronary ischemia
- Arterial occlusion/thrombus, near the puncture site
- Arterial occlusion/thrombus, remote from puncture site
- Arterial occlusion/restenosis of the treated vessel
- Arteriovenous fistula
- Arrhythmia
- By-pass Surgery
- Death related to procedure
- Death unrelated to procedure
- Drug reactions to antiplatelet drugs/contrast media
- Embolization, arterial
- Embolization, stent
- Fever
- Hematoma bleed, remote site
- Hematoma bleed, puncture site
- Hypotension/hypertension
- Incompatibility of the local anesthetics
- Intimal injury/dissection
- Ischemia/infarction of tissue/organ
- Local infection

- Malposition (failure to deliver the stent to the intended site)
- Pulmonary embolism
- Perforation or Dissection of the vessel
- Pseudo aneurysm
- Renal failure
- Radiopaque material induced nephropathy
- Restenosis
- Rupture of the retroperitoneum or an adjacent organ
- Septicemia/bacteremia
- Stroke/cerebrovascular events
- Vascular Complications (hematoma, arterial fistula, aneurysms spuria and arterial thromboses)
- Vasospasm
- Venous occlusion/thrombosis, remote from puncture site
- Venous occlusion/thrombosis, near the puncture site
- Sepsis/infection
- Stent migration
- Stent fracture
- Stent misplacement

## IX. DIRECTIONS FOR USE

### Pre-Deployment Procedure

1. Inject Contrast Media  
Perform an angiogram using standard technique.
2. Evaluate and Mark Target Site  
Fluoroscopically evaluate and mark the target site, observing the most distal diseased or obstructed segment.
3. Select Stent Size  
Measure the length of the target lesion or stricture to identify the appropriate length of stent required. Ensure that the stent is long enough to permit the area proximal and distal of the lesion or stricture to be covered by the stent.  
Identify the diameter of the reference vessel (proximal and distal to the lesion or stricture). To ensure secure placement, refer to the stent size selection table for proper sizing scheme. Refer to product labeling for stent length.

Reference Vessel Diameter	Unconstrained Stent Inner Diameter
4.0 - 5.0 mm	6.0 mm
5.0 - 6.0 mm	7.0 mm
6.0 - 7.0 mm	8.0 mm
7.0 - 8.0 mm	9.0 mm
8.0 - 9.0 mm	10.0 mm
9.0 - 10.0 mm	11.0 mm

Table 1

4. Prepare the Peripheral Vascular Self Expanding Stent System
  - a) Open the box and remove the pouch containing the stent system.
  - b) Carefully inspect the pouch for damage to the sterile barrier. Then, peel the pouch open and remove the tray containing the stent system. Extract the stent system from the tray and check the following:
    - i) Verify that the shipping lock (Figure 3) is still secure in the stent system handle.
    - ii) Examine the stent system for any damage. If it is suspected that the sterility or performance of the device has been compromised, the device should not be used.
  - c) Visually inspect the distal end of the stent system to ensure that the stent is contained within the sheath. DO NOT use if the stent is partially deployed.
  - d) Flush the inner lumen of the device with saline prior to use.
  - e) Wipe the usable length portion of the stent system with gauze soaked with saline.

### Stent Deployment Procedure

1. Insert Introducer Sheath and Guidewire
  - a) Gain access at the appropriate site utilizing a 6F (2.0 mm) (or larger) introducer sheath.
  - b) Insert a guidewire of appropriate length (see table) and diameter across the stricture to be stented via the introducer sheath.

Usable Catheter Length of delivery device	Recommended guide wire length
120 cm	300 cm
80 cm	260 cm

Table 2

2. Dilate Stricture  
If the physician deems that predilation is required, standard techniques (e.g. PTA balloon catheter) may be used. While maintaining site access with a guide wire, remove the PTA balloon catheter from the patient.  
**Caution:** During dilation, DO NOT expand the PTA balloon such that dissection complication or perforation could occur.
3. Introduce the Peripheral Vascular Self Expanding Stent System
  - a) Advance the device over the guidewire through the sheath introducer.  
**Note:** If resistance is met during stent system introduction, the system should be withdrawn and another system should be used.  
**Caution:** Always use an introducer sheath for the implant procedure to protect the vasculature and the puncture site. A 6F (2.0 mm) (or larger) introducer sheath is recommended.
  - b) Position the tip of the delivery device catheter past the target site.
  - c) Pull back the delivery device until the distal and proximal ends of the stent are in position so that they are distal and proximal to the target site.

Figure 4 – Slack in stent system

- d) Remove slack (see marked area in Figure 4) from the stent system held outside the patient.

**Caution:** Any slack in the stent system (outside the patient) could result in deploying the stent beyond the target site.

4. Deploy Stent
  - a) Verify that the distal and proximal stent ends are distal and proximal to the target site.
  - b) Confirm that the introducer sheath is secure and will not move during deployment.
  - c) Remove the shipping lock (Figure 2).
  - d) Hold the handle in one hand by grasping the device by placing the handle in the palm of your hand in a pistol like fashion.
  - e) Leave your thumb bent and in a bent fixed position as in Figure 5.
  - f) It is important to maintain this position throughout the deployment of the device as a modification of this position will enhance the chance of operating the device in a syringe like fashion which may cause the deployment of the prosthesis to be more distal than optimal.
  - g) Please note that a modification of your thumb to a position of that in Figure 6 is contra indicated and should not be used.
  - h) With your hand firmly around the handle assembly, and your thumb in a cocked position, your index finger on the top of the first handle assembly as show in Figure 7 above and your middle and ring finger firmly resting on the middle assembly as shown in Figure 8 slowly pull back on the 1<sup>st</sup> trigger assembly with a firm non-jerking motion.
  - i) As the prosthesis begins to deploy you will have the tactile sensation of the outer sheath of the device slowing retracting to release the prosthesis.
  - j) Verify the position of the device using fluoroscopy.
  - k) Continue retracting the outer tube by slowly pulling the first trigger of the handle.
  - l) For longer prosthesis a second trigger assembly should be engaged in the same fashion as the first.
  - m) Once your index, middle, and ring finger are firmly on the handle assembly begin to deploy the second trigger.
  - n) Pull the second trigger of the handle until the stent is released completely.
  - o) It is contraindicated to push the hand assembly or to try to reconstrain the stent while using this device and any attempt to do so may cause the stent to be placed in a less than optimal location.


Figure 5


Figure 6


Figure 7


Figure 8

5. Post Stent Placement
  - a) Remove the delivery device from the body.
  - b) If additional stent-to-vessel apposition is desired, select a balloon catheter that matches the size of the reference vessel, but that is not larger than the stent diameter itself.
  - c) Remove the guidewire and introducer sheath from the body.
  - d) Close entry wound as appropriate.
  - e) Discard the delivery device, guide wire, and introducer sheath.

**Note:** Physician experience and discretion will determine the appropriate drug regimen for each patient.

## X. MAGNETIC RESONANCE IMAGING (MRI COMPATIBILITY)

Non-clinical testing has demonstrated that the Peripheral Vascular Self Expanding Stent System is MRI-compatible. Patients with the Peripheral Vascular Self Expanding Stent System can be scanned safely, immediately after placement of this implant, under the following conditions:

- Static magnetic field of 3-Tesla or less
- Spatial gradient field of 720 Gauss/cm or less
- Maximum specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning.

In non-clinical testing, the Peripheral Vascular Self Expanding Stent System produced a temperature rise of less than or equal to 1.4°C at a maximum specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning in a 3-Tesla MR system (Excite, Software G3,052B, General Electric Healthcare, Milwaukee, WI).

MR image quality may be compromised to the area of interest is in the exact same area or relatively close to the position of the Peripheral Vascular Self Expanding Stent System. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

## XI. REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

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